


The XXIXth Annual Occupational Ergonomics and Safety Conference
Seattle, Washington, USA
June 1-2, 2017

*The Sling Specialist
In Safe Patient
Mobilization*



SLING CATALOG
Reusable and Wipeable Series
January 2017

800.273.5749
www.alphamodalities.com
info@alphamodalities.com



ISBN: 97819384965-5-4



Color Coding

Guide to Color Coding for Mobility Assessment

RED Max Assist; Unable or does not cooperate

Page 1.....A-TSL Series
Page 2.....A-Seat Series
Page 3.....A-TTurner Series
Page 4.....A-Pannus/A-PannusSW

ORANGE Minimal to Moderate Assist; Stand by Assist
or Contact Guard Assist

Page 5.....A-SitStandSW
Page 6.....A-WalkVSW
Page 7.....A-BariV/A-BariVSW
Page 8.....A-Limb
Page 9.....A-WedgeSW

YELLOW Minimal to Moderate Assist; Stand by Assist

Page 10.....A-BedLadderSW

A-TSL Series

Intended use: Turns, boosting up in bed, lateral transfers, and recumbent transfers.



Remains under patient



Soft to the touch

- Washable
- Latex Free

Item/Color	Size	SWL
A-TSLITF8-32	32"W 60"L/ 81.3cm x 152.4cm	660 Lb/ 300 Kg
A-TSLITF8-45	45"W 86"L/ 114.3cm x 218.4cm	1000 Lb/ 454 Kg
A-TSLITF8-45/100	45"W 100"L/ 114.3cm x 254cm	1000 Lb/ 454 Kg
A-TSLITF8-50	50"W 86"L/ 127cm x 218.4cm	1000 Lb/ 454 Kg
A-TSLITF8-60	60"W 86"L/ 152.4cm x 218.4cm	1000 Lb/ 454 Kg
A-TSLMEF8-45	45"W 86"L/ 114.3cm x 218.4cm	1000 Lb/ 454 Kg
A-TSLMEF8-50	50"W 86"L/ 127cm x 218.4cm	1000 Lb/ 454 Kg
A-TSLMEF8-60	60"W 86"L/ 152.4cm x 218.4cm	1000 Lb/ 454 Kg

IT- proprietary new innovation in sling fabric, incorporating three key features: breathability, wicking and antimicrobial
 ME-open weave mesh, breathable, multiple colors



A-Seat Series

Intended use: Seated transfers to and from bed, chair, toilet and wheelchair with head support.

Washable series-
piping color indicates
size



Option to add
amputee strap
system

Mesh option
for
bathing

Wipeable series-
outside shell color
indicates size



Two dignity
straps on
wipeable

- Washable and wipeable* option
- Latex Free

Size/ Trim Color	Recommended Weight range
A-Seat-S	31-75 Lb / 14-34 Kg
A-Seat-M	75-145 Lb / 34-66 Kg
A-Seat-L	145-220 Lb / 66-100 Kg
A-Seat-X	220-440 Lb / 100 -200 Kg
A-Seat-XX	440-800 Lb / 200-363 Kg

*Wipeable fabric is non breathable
Wipeable fabric does not come in Size XX

A-TTurner Series

Intended use: To turn patient and to hold patient onto their side.



Fits into lower
back of patient

Bariatric size
available



- Washable and wipeable* option
- Latex Free

Item	SWL
A-TTurner	330 Lb/ 150 Kg
A-TTurnerB	660 Lb/ 300 Kg
A-TTurnerBSW	660 Lb/ 300 Kg

*Wipeable fabric is non breathable

A-Pannus/A-PannusSW

Intended use: To assist in moving patient's panniculus away from the body for cleaning skin, wound checks, or dressing changes.



- Washable and wipeable* option
- Latex Free

Item	SWL
A-Pannus	1000 Lb/ 454 Kg
A-PannusSW	1000 Lb/ 454 Kg

*Wipeable fabric is non breathable

A-SitStandSW

Intended use: Used with a powered mobile lift unit to assist a patient to a standing position.

Inner safety belt
that is adjustable and
nonslip



Color coded
tab indicates
size



Clip attachment
options



- Wipeable*
- Latex Free

Item	Waist	SWL
A-SitStandSW-S	28"-44"/ 71.1-111.8cm	660 Lb/ 300 Kg
A-SitStandSW-M	38-55"/ 96.5-139.7cm	
A-SitStandSW-L	45-64"/ 114-162.6cm	
A-SitStandSW-XL	55-73"/ 139.7-185.4cm	

*Wipeable fabric is non breathable

A-WalkVSWWL

Intended use: Ambulation and gait training



Groin straps
unclip for toileting
option

Sling color
indicates size



- Wipeable*
- Latex Free

Item	Chest Size	SWL
A-WalkVSWWL-Peds	20"-27" / 50.8--68.58cm	660 Lb/ 300 Kg
A-WalkVSWWL-S	32"-39" / 81.3--99.1cm	
A-WalkVSWWL-M	40"-44" / 101.6-111.8cm	
A-WalkVSWWL-L	42"-47" / 106.7-119.4cm	
A-WalkVSWWL-XL	49"-54" / 124.5-137.2cm	

*Wipeable fabric is non breathable

A-BariV/A-BariVSWWL

Intended use: Ambulation and gait training



- Washable and wipeable* option
- Latex Free

Item	Waist	Groin	SWL
A-BariV/ VSWWL-100	20-40"/51-102cm	12-20"/30-15cm	220 Lb/ 100 Kg
A-BariV/ VSWWL-200	40-70"/102-178cm	40-59"/102-150cm	440 Lb/ 200 Kg
A-BariV/ VSWWL-300	63-94"/160-239cm	63-88"/160-226cm	660 Lb/ 300 Kg
A-BariV/ VSWWL-400	74-112"/188-284cm	74-102"/188-259cm	880 Lb/ 400 Kg

*Wipeable fabric is non breathable

A-Limb

Intended use: To support limbs for dressing, wound, and/or back care.



Clip or loop attachment option

Custom sizing available



- Washable and wipeable* option
- Latex Free

Item	Size	SWL
A-LimbF8	33"x8" / 83.8cm x 20.3cm	1000 Lb/ 454Kg

*Wipeable fabric is non breathable

A-WedgeSW

Intended use: Positioning wedge to keep patient on their side; prevent sliding down in bed.

Two size options



Two inner materials options:
Foam (FM)-Firm
Bead (BN)-Pliable

- Wipeable*
- Latex Free

Item	Size
A-WedgeSWFM-M	16"L (40.64cm) Positioning
A-WedgeSWFM-L	23"L (58.42cm) Positioning
A-WedgeSWBN-M	16"L (40.64cm) Contour
A-WedgeSWBN-L	23"L (58.42cm) Contour

*Wipeable fabric is non breathable

A-BedLadderSW

Intended use: Allow patient to sit themselves up and lift limbs with minimal to no assistance.



Soft, pliable grips



MRI
Compatible



- Wipeable*
- Latex Free

Item	Size	SWL
A-BedLadderSW	63" / 160.02cm	600Lb/ 272Kg

*Wipeable fabric is non breathable

Compatibility Statement

All lift motors perform a basic function, which is to go up and down. The interface between the lift, hanger bar and sling is critical for safety and patient comfort.

We test our slings with various manufacturers' hanger bars for compatibility to ensure that they are safe to use and are comfortable for the patient.

We guarantee compatilbly with the following manufacturers:

Guldmann AS
Arjo Huntleigh
Prism Medical
Handicare
Tollos
EZ Way
Hill Rom Liko
Etac-Molift



Customizations

Sling Fabric Colors




Also available
in white or
other custom
colors



Sling Label Options

Multiple attachment options



Cleaning and Warranty



Cleaning instructions for Wipeable products:

- Clean with approved facility disinfectant for "Non Critical Devices", as per FDA guidelines
- Do not soak in solution of chlorine based bleach solution
- Do not autoclave or steam clean



Cleaning instructions for Washable products:

- Wash Temp not to exceed 185F / 85 C
- Drying- Medium Heat not to exceed 165F / 74C
- Oxygen Based Bleach Recommended
- Max Load 850lbs Cycle

Warranty

Within the first Three (3) Years of use or production date, we will repair or replace, at no charge, any sling found to be defective in manufacturing or damaged in use /cleaning. Sling must be returned to us clean, freight prepaid and attached to Return Authorization Number.

ALL PRODUCTS ARE LATEX FREE



ISBN: 97819384965-5-4

3110 Judson Street Ste 149, Gig Harbor WA 98335 • 800.273.5749 • www.alphamodalities.com

